

Our product range



Safety Lancet



Lancing Device



Universal Lancet



Blood Lancet

NEW: Neonate Heel Incision Devices



Steri-Heel™
Premie



Steri-Heel™
Newborn

VitalCare UK

Jubilee House
The Drive, Great Warley
Brentwood, Essex
CM13 3FR, UK

t: +44 (0) 1277 725188

f: +44 (0) 1277 725001

e: info@vitalcare.co.uk

w: vitalcare.co.uk


VitalCare

Auto Safety Lancets



Let's get to the point

High quality lancets at great prices.



Better value for your money

The *AUTO* is a small, disposable safety device which administers a prick to the finger, to obtain a drop of blood.

The *AUTO* is by far the most cost-efficient safety lancet available on the NHS contract. You could make very substantial savings simply by switching to VitalCare.

The *AUTO* is a worldwide success due to its simple yet effective patented trigger technology which delivers a safe, virtually pain-free finger prick – every time.

The highly polished, finely ground needle performs the prick in a split second, before retracting back into the device – for good. The *AUTO* is entirely self-disabling, so there is no risk of needle-stick injury or cross-infection.

We can cater for all types of skin and patient groups.

Model (Auto)	N P C	Needle size
Junior pink	FTM065	26g x 1.8mm
Soft blue	FTM066	21g x 1.8mm
Normal beige	FTM067	21g x 2.4mm
Extra green	FTM069	21g x 3.0mm

Model (Baby)	N P C	Blade size
Preemie orange	FTM254	0.85mm (D) x 1.75mm (L)
Newborn yellow	FTM253	1.0mm (D) x 2.5mm (L)

Why pay substantially more for the same outcome?

The price table below shows just how much your Trust could be saving:

Brand	per 100 devices
Vitalcare AUTO	£6.48*
Sarstedt	£7.54**
Surgilance	£8.08*
Unistik 3	£7.65*
Haemolance	£8.89**
Safe-T-Pro Plus	£9.54
Medlance II	£10.00
Glucocet / Minilet	£11.84

Prices taken from NHSSC Blood Lancets and Associated Products – Product listing Summary Oct 2011.

* Band 1 price per box of 100pcs

** Band 1 price per 100pcs (equivalent)

